

DATA SHEET

VETRAP 3M

DOSAGE AND METHOD OF USE

Under no circumstances should the band be applied with a higher tension than 50% of the maximum tension.

SPECIAL WARNINGS

- **Vetrap 5 cm:** blue, red, green
- **Vetrap 7,5 cm:** blue, red, green
- **Vetrap 10 cm:** blue, fuchsia, yellow, black, green, red



PACKAGING

- **Vetrap 5 cm:** box of 200 pieces (can be sold individually)
- **Vetrap 7,5 cm:** box of 150 pieces (can be sold individually)
- **Vetrap 10 cm:** box of 100 pieces (can be sold individually)

DESCRIPTION

Elastic bandage for restraint and maintenance. Lightweight, strong and porous, suitable for wound protection, holding dressings in place and tail dressing.

- Adheres to itself without the need for additional fixings
- maintains mobility
- easy to remove, does not stick to the hair
- allows the bandaged area to breathe
- water- and sweat-proof
- if sufficient pressure is applied, it can be used to stop bleeding while remaining flexible enough not to impair blood circulation

