### **DATA SHEET**

# **MASH PROTECT**



# **COMPOSITION**

Oats • maize flakes • dehulled oats •

linsed meal • sunflower seed meal • alfalfa meal • wheat flakes • dextrose • yeasts • whey powder • chicory inulin • yeast cell walls from Saccharomyces cerevisiae • ground Psyllium husk • product from Lactobacillus spp rich in proteins

#### ADDITIVES/KG

Digestibility enhancer: Saccharomyces cerevisiae CNCM I-4407 (4b1702) 1,3X1010 CFU

## DOSAGE AND METHOD OF USE

Rehydrate 600 g (1 sachet) of Mash Protect with warm water (do not use more than 1 litre).

- Cover the container to allow the mash to swell
- · Allow to cool before administering

Mash Protect can be routinely administered once a week as a meal replacement

up to a maximum of 3 times a day (e.g. after surgery)

# **PACKAGING**

Box of 9 sachets (600 g each)

# **DESCRIPTION**

Ready-to-use complete feed, packaged in practical single doses.

Mash Protect contains high quality raw materials consisting of selected cereals to which are added ingredients useful for maintaining intestinal habitat and the functioning of digestive processes.

# **PRODUCT CODE**

EAN13: 8055320270222

